


NIO™ Infant

Manual, quick, safe, easy vascular access from 2.3kg to 3 years old

- Facilitates rapid vascular access in infant patients
- Requires no batteries or maintenance
- Pocket-sized and lightweight
- Fixation dressing helps to prevent dislodgement
- Sterile, single-use device
- Stepped Needle®*
- Safety needle cap
- 5-year shelf life



Vascular Access in a Flash

The NIO™ Infant provides safe, rapid intraosseous access in the proximal tibia as an alternative to IV access in emergency medical situations. The NIO Infant is a manual IO device with an innovative Stepped Needle® safety design which provides graduated penetration into the bone marrow cavity to help prevent overpenetration.

The NIO Infant is ready to use out of its packaging. The device requires no batteries, no assembly, and no extra parts, allowing for rapid vascular access to be achieved.

*Patent pending

NIO™ Infant

Product Highlights

Pocket-sized, lightweight device facilitates rapid vascular access

Safety needle cap decreases the risk of accidental needle injury

Stepped Needle prevents overpenetration

Fixation dressing helps to prevent dislodgement



Specifically designed for infants—meet the NIO Infant:

- **Lightweight, ergonomic handle** ensures simple, safe manual insertion of the needle into the bone marrow cavity
- Stepped Needle design allows the user to **manually control penetration** into the bone
- The Stepped Needle provides sensory feedback in the form of resistance to the user, indicating that the needle is properly positioned in the bone marrow cavity, thus, **preventing overpenetration**
- Fixation dressing features adaptable clasp which accommodates various penetration depths to **prevent dislodgement**

Specifications

| | |
|---------------------------------------|--|
| Weight | 56g |
| Packaging | Packed individually in a peel back blister pack with fixation dressing |
| Safety Mechanism | Stepped Needle® 14G to 18G |
| Needle Gauge | 18G |
| Needle Effective Length | 15mm (.6in) |
| Needle Penetration Depth | The cannula cortex penetration section (18G) is 2.6mm (.1in) |
| Shelf Life* | 5 years* |

Specifications are subject to change without notice.
*From date of manufacture

PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com

www.persysmedical.com

