# MEDIEJECT II

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

FIG 3

### 1. FOREWORD

GCE medical suction ejectors are medical devices classified as class IIa according to the Medical Device Directive 93/42/EEC.

Their Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon standards EN ISO 10079-3 and EN 1789, as amended.

### 2. INTENDED USE

MediEject II is a medical device designed to generate continuous vacuum with continuous regulation. It is an ejector that works on the principle of Venturi effect. The propellant is compressed air or oxygen. The device is connected in a system with a collecting container or with multiple containers and a filter connected by suction hoses.

The system may be fixed or portable for use inside healthcare facilities and also in an ambulance.

The product is not intended for use during thoracic drainage.

The suitability for a specific procedure shall be decided by appropriately qualified medical staff.

### 3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

Keep the product and its associated equipment away from:

- All sources of heat
- Flammable materials
- Oil or grease (including all hand creams)
- Water
- Dust.

The product and its associated equipment must be prevented from tipping over and falling.

Always maintain oxygen cleanliness standards.

Use only the product and its associated equipment in a well ventilated area.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-10/+40 °C		-20/+60 °C
	10/70 %		10/100 %
	950/1100 mbar		600/1200 mbar

Before initial use the product shall be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdraw from operation (for transport and storage). Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

## 4. PERSONNEL INSTRUCTIONS

The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

**⚠ Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.**

## 5. PRODUCT AND FUNCTION DESCRIPTION

MediEject II suction ejector includes 3 basic variants:

- MediEject II with probe\*
- MediEject II T-slide
  - with NIST
  - with HOSE and probe\*
- MediEject II Rail clamp\*
  - with NIST
  - with HOSE and probe\*

\* according to national standards

It mainly consists of a probe connector (gas source side) (A), an integrated push-pull on/off valve and control handwheel (B), a vacuum indicator (C), a hose nipple (patient side) (D) and silencer (E).

### A. Probe Connector (medical pipeline or mobile gas source side)

The product is fitted to the medical gas pipeline system or mobile gas source by a oxygen/air type male probe. The probe is directly mounted on the ejector or on the hose via rail mounting system.

Hose may be connected via NIST connector and does not have to be part of the delivery, depending on the ordered variant.

Refer to appendix Nr 2 to get information about connection and release method.

#### B. Integrated handwheel

MediEject II is equipped with an integrated controller with a push-pull on/off valve for closing gas source and vacuum adjustment by means of a knob.

The suction ejector is on when the controller is pulled up and off when it is pressed down. A green label indicates that the controller is on (ON). The integrated controller allows the user to set vacuum at the patient side in the range from 0 to 80 kPa. The vacuum level at the patient side can be increased by rotating the control knob counter clockwise. The vacuum level at the patient side can be decreased by rotating the control knob clockwise.

#### C. Vacuum indicator

MediEject II is equipped with a vacuum gauge (manometer) indicating the value of vacuum on the patient side during flow.

 **If the patient side is without flow, the vacuum value will remain at maximum when the knob is turned clockwise. See non-return valve description.**

#### D. Hose nipple (patient side)

The suction ejector has a hose nipple at the patient side that enables to connect the suction ejector with associated equipment such as the suction hose.

#### E. Silencer

The silencer is intended to reduce a noise produced by the device when creating vacuum.

#### F. Non-return valve

The non-return valve is located in the ejector before the hose nipple. It is a safety feature that prevents pressure from entering the patient side in the event of a malfunction (such as a clogged silencer). Likewise, the non-return valve prevents vacuum drop at the hose nipple, if the vacuum is lowered by rotating the control knob clockwise in the absence of flow. This drop is not indicated on the vacuum indicator. See vacuum indicator.

## 5.1. TECHNICAL DESCRIPTION

Power gas:	O <sub>2</sub> , AIR *
Inlet pressure:	400 - 500 kPa
Max. vacuum:	20-80 kPa <a href="#">1)</a> <a href="#">2)</a> <a href="#">3)</a>
Max gas consumption at inlet pressure 400 kPa:	25 l/min

Free flow suction at inlet pressure 400 kPa:	25 l/min <sup>4)</sup>
Noise - close/open suction:	35/45 dB
Outlet connection:	Hose nipple - inner diameter 6 mm

EN

\* For gas used, see Appendix Nr. 1

The parameters for the respective variant can be found in Appendix No. 1.

According to EN ISO 10079-3 it is a device with High vacuum level at least 60 kPa below atmospheric pressure and with High flowrate at least 20l/min.

1) Reference conditions: 20°C;1013hPa

2) The achievable vacuum depends on atmospheric pressure.

Every 100m increase in altitude causes a decrease in the achievable vacuum of about 1.4 kPa

3) Total accuracy of vacuum level is  $\pm 8$ kPa


4) The control range of free suction flow is 10-25l/min

## 6. OPERATIONS

### 6.1. CONNECTION & USE OF INLET (PATIENT SIDE)

MediEject II is used with the optional associated equipment specified in Chapter 8.

**⚠ MediEject II installation must be in accordance with EN ISO 10079-3 and EN 1789 as amended. For connection, see the figure below.**

1. Suction ejector

2. Filter

3. Suction hose

4. Collection container\*

**⚠ During suction, MediEject II must be used with one or more collecting containers, which must include an overflow protection device to prevent liquid from entering the suction ejector.**

-  The air leaving the collection container must pass through a microbiological filter before entering the suction device.
-  The associated equipment or medical devices connected to the suction ejector must be fully compatible with the characteristics and performance data of the MediEject II suction ejector.

\* In some markets, it is recommended to use a safety container installed between the vacuum source and the collection container before the ejector.

## 6.2. BEFORE USE

### Visual Inspection before fitting

- Check if there is visible external damage on the suction ejector and its associated equipment (including product labels and marking). If it shows sign of external damage, remove the products from service and identify its status.
- Visually check if the ejector and its associated equipment are contaminated; if needed, apply the cleaning procedure according to cleaning procedure here under.
- Check if the service or the disposal time has not been exceeded, using GCE or owner's date coding system. If service or disposal time has been exceeded, remove the product from service & suitably identify its status.

**Fitting to medical gas pipeline system - terminal units/mobile gas source**  
Refer to Appendix Nr. 2 to get information about connecting / disconnecting procedure.

For ejector fitted with a rail clamp and a flexible hose, clip the suction ejector on the rail prior to connect the quick coupler.

### Functional test before use and Leakage test

- After connecting the medical pipe line system terminal unit/portable gas source and associated equipment, turn the vacuum control knob to maximum position (counter clockwise).
- Set the suction ejector ON by pulling up the integrated handwheel.
- Check that there is a significant vacuum level at patient side inlet by locking the patient side inlet with your finger.
- Set the suction ejector OFF by pushing the integrated handwheel and check if vacuum level is without change at locked patient side inlet by your finger.
- Set vacuum level to minimum by turning the integrated handwheel it clockwise.

-  Before connecting any associated equipment to the product, make sure that the push-pull on/off valve is in the closed position and that the patient is not connected.

-  If any leaks are found, check the connection and repeat the leak test.

- ⚠** If a leakage is detected again, apply the “After use” procedure below and return the valve for service or repair.

### 6.3. OPERATIONS

- Set the suction ejector ON by pulling up the integrated handwheel.
- Lock the patient side inlet with your finger.
- Set the vacuum level to the needed level for your therapy by turning the integrated handwheel in the anti-clockwise direction and checking the vacuum indicator (gauge).

- ⚠** If vacuum value is reduced in the absence of flow, the vacuum indicator remains at maximum. See the description of the non-return valve.

- The suction ejector is ready to be used.

### 6.4. AFTER USE

- Before disconnecting ejector from the medical pipe line system terminal unit/mobile gas source:
  - turn it OFF, by pushing down the integrated handwheel and set vacuum level to minimum by turning the integrated handwheel clockwise
  - remove suction tubing connection from patient side

- ⚠** Follow the associated equipment manual. For example: If the collecting container is overfilled, extracted fluids are sprayed on / into the patient when the vacuum source is abruptly interrupted.

## 7. ACCESSORIES

Not available.

For closer information about the possible optional equipment refer to the section 8.2 associated equipment.

## 8. CLEANING

### 8.1. SUCTION EJECTOR


- ⚠** MediEject II must be cleaned and disinfected externally after every patient.

If there is risk that suction ejector is exposed to contamination following cleaning shall be used:

Remove general contamination with a soft cloth damped in oil free oxygen compatible soap water and clean with deionised water.

Disinfection can be done with an alcohol-based preparation (moistened wipe). If another cleaning agent is used, check that it is not abrasive and is compatible with the product materials (including labels) and gas (a suitable cleaning agent - e.g. Meliseptol).

- ⚠ Do not use cleaning solutions containing ammonia!
- ⚠ Do not immerse the product in any liquid.
- ⚠ Do not expose to high temperature (such as autoclave).
- ⚠ Clean immediately after use.
- ⚠ Clean in position as shown in this chapter. Do not turn upside down.
- ⚠ Do not use spray to apply the cleaner, as spray can penetrate the inside of the product and cause contamination or damage.
- ⚠ Do not use pressure water as it may damage or contaminate the product.
- ⚠ Do not use the product under any circumstances if the internal parts of the product have been contaminated. It must be decommissioned.

There is a ventilation opening for the vacuum indicator located below the face plate at the bottom of the vacuum indicator. See the mark „-“.

- ⚠ Keep the ventilation area clean and dry.

## 8.2. ASSOCIATED EQUIPMENT

The used associated equipment shall comply with EN ISO 10079-3 and EN 1789 as amended and with appropriate Instruction for use:

- Microbiological filter: Use a filter with the required connection.
- Suction hose: Use suction hoses with an inner diameter of 6mm.
- Collection container: Use the collection container with the required connection.
- Safety container: Use safety containers with the required connection

- ⚠ Before connecting any associated equipment or medical device to the MediEject II, always check that they are fully compatible with connection features & performances of the product.

## 9. MAINTENANCE

### 9.1. LIFE TIME

#### Life time and waste management

Maximum life time of the product is 10 years.

At the end of the product's life time (10 years maximum), the product must be withdrawn from service. The provider of the device shall prevent the reuse of the product and handle the product in compliance with "Directive of European Parliament and Council 2008/98/EC on waste".

In accordance to Article 33 of REACH GCE, s.r.o. as responsible manufacturer shall inform all customers if materials containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC). The most commonly used brass alloys used for bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-468-6, CAS no. 7439-92-1. The lead will not be released to the gas or surrounding environment during normal use. After end of life the product shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact to environment and health.

To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE product.

#### **Functional test**

Check the function of the device before every use (chapter 6.2).

#### **Serial number and date of production**

Form of nine digit serial number stamped on the product is following:

YY MM XXXXX

YY: year of production

MM: month of production

XXXXX: sequence number

Example: serial number 090300521 shows the product produced in March 2009, with sequence number 521.

## **9.2. REPAIRS**

The repairs shall be carried out by a GCE authorised person only.

Any product sent back to a GCE authorised person for repair shall be properly packaged. The purpose of the repair has to be clearly specified. A short description of fault and any reference to a claim number might be helpful.

Some repair activities concerning to the replacement of the damaged or missing components can be carried out by the owner of the product. The following parts can be replaced only:

- Labels <sup>1)</sup>
- Silencer (E)
- Gauge protective glass (C)


1) Before stick Front label clean area under gauge from old glue and clean groove if the product is equipped.

 **Use only original GCE components.**

 **All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.**

 **GCE MediEject II valve do not contain any other provider replaceable parts.**

## 10. GLOSSARY

	Consult instruction for use		Suitable for Hospital care use
	Caution		Suitable for Emergency care use
	Keep away from heat and flammable material		Serial number
	Keep away from oil and grease		Catalogue number
	Humidity limit		Batch code
	Temperature limit		Fragile, handle with care
	Keep dry		Manufacturer
	Date of manufacture		Outlet parameter
	Use by date		Inlet parameter
	High vacuum/ High flow		

EN

## 11. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product). The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

**APPENDIX:**

Nr. 1: Technical description and performance

Nr. 2: Quick coupling feature and connecting/disconnecting procedure

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© GCE, s.r.o.