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Dafilon®

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Dafilon®

Description

Dafilon® is a sterile, monofilament, non-absorbable surgical suture material produced from polyamide. Dafilon® is available in blue (polyamide 6/66), coloured with Copper Phthalocyanine (C.I. 74 160) or in black (polyamide 6/6), coloured with Hematin Logwood Black (C.I. 75290) to make it more readily recognizable, but is also available undyed in the natural beige colour. Dafilon® fulfils all the requirements of the European Pharm. and the United States Pharm. current edition for sterile, non absorbable sutures.

Indications

Dafilon USP sizes from 5 to 6/0 are intended to be used for soft tissue approximation in skin closure and oral surgery. Dafilon USP sizes from 8/0 to 11/0 are intended to be used for

for soft tissue approximation in Microsurgery, Ophthalmology, Microsurgical repair and Neurosurgery.

Dafilon is intended to be used for adult and paediatric patients who do not meet any contraindication.

Mode of action

Suture materials are used primarily for adaptation of the wound edges to render possible an undisturbed wound healing. When Dafilon® suture materials are employed, a minimal inflammatory tissue reaction will occur followed by a gradual encapsulation of the suture material by fibrous connective tissue. Like all polyamide materials, over a long period of time Dafilon® in vivo is subject to a gradual loss of the total suture strength through the action of body enzymes.

Contra-indications

On account of the gradual loss of strength, which can occur in vivo over a longer period, the use of Dafilon® is contraindicated when a permanent retention of suture strength is of importance.

Dafilon USP sizes from 5 to 6/0 are contraindicated for cardiovascular surgery and neurosurgery.

It is contraindicated to patients with hypersensitivity or allergy to the suture material.

Mode of Application

Dafilon should be used applying the standard surgical suturing and knotting techniques, taking into account the surgeon's experience with the respective surgical procedure. Sutures should be selected depending on the tissue to be sutured, the specific duration of wound support needed, the size of the wound, the patient's condition and the specific suturing technique. When Dafilon is used for skin closure and oral surgery, and being a non- absorbable material, it is strongly recommended to remove the suture once the healing period is achieved, and no later than 30 days.

Warning

Dafilon® must not be re sterilized. Opened unused or damaged packs shall be discarded.

Dafilon must not be used: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

Precaution

The user should be familiar with surgical suturing techniques, before employing Dafilon® suture materials. When working with suture materials great care should be taken to ensure that the use of surgical instruments, such as tweezers and needle holders, does not cause the material to be damaged by being pinched or kinked. Care should be taken to avoid damaging the needle when using the suture material. Always grasp the suture in a section 1/3 to 1/2 of the distance from the fiber attachment end to the needle point, never at the end where the fiber is attached or the needle point. Grasping the needle at the area of its point could impair the penetration performance and cause a fracture of the needle. Grasping the needle close to the fibre attachment end could cause bending, breakage. Reshaping needles should be avoided and may result in a loss of their strength and resistance towards bending and breaking.

Exertion of inappropriate tension on the needle-fiber fixation can lead to a detachment of the fiber from the needle. Care must also be taken to avoid needle stick injuries when handling with surgical needles. Discard the needles after use in specially intended containers. Care should be taken that the knots are positioned properly and adequate knot security is given.

Side effects

As for every other suture material, prolonged contact with body fluids such as urine and bile, can lead to lithiasis. The following side effects may be associated with the use of this product: tissue injury, transitory local irritation, transient inflammatory foreign body reaction, enhanced bacterial infectivity, wound dehiscence, hypertrophic scar/granuloma formation, pain, seroma, hematoma, edema, tract bleeding, increased risk of aneurism and stenosis.

In rare cases hypersensitivity or allergic reaction to the suture material have been observed.

Sterilization

Dafilon sutures can be sterilized by ethylene oxide gas or gamma radiation.

Storage

No special storage conditions are required.

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