Veroval® Netzgerät

Deutsch

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CH IVF HARTMANN AG
8212 Neuhausen
www.veroval.ch

Made under the control of PAUL HARTMANN AG by
Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road
European Industrial Zone, Xiaolan Town
528415 Zhongshan City, Guangdong Province
Made in P.R.C.

EC REP Donawa Lifescience Consulting Srl
Piazza Albania, 10
00153 Rome
Italy

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Veroval® AC/DC Adaptor

English


Operating instructions
AC/DC adaptor for HARTMANN Veroval® Blood Pressure Monitors


Product description
The main part plus converts the alternative voltage of 100-240 V to a safety extra low voltage. The device is tested according to the medical standard IEC 60601-1 and is intended for use with compatible Veroval® blood pressure monitors.

- **Safety and precautionary measures**
- Read the operating instructions before initial use.
- Observe the Veroval® blood pressure monitor instructions for use.
- Only connect the power supply unit to power sources that are suitable for this purpose. Ensure the polarity is correct.
- The power supply unit must never be exposed to rain, moisture or high heat.
- Only operate power supply unit indoors.
- Do not use the power supply unit when damaged.
- Do not open or modify the power supply unit. Repairs may only be carried out by authorised technical personnel.
- Do not cover the power supply unit. Ensure sufficient cooling. Do not expose to direct sunlight.
- Keep power supply unit and cables away from flammable items.
- Only clean with a damp cloth. Do not use any cleaning products. Do not sterilise device.

- **Notes on electromagnetic compatibility**
- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- It may only be possible to use the device to a limited extent in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly beside other devices or stacked on top of other devices, as this could result in faulty operation. However, if it is necessary to use the device in the manner described above, this device and the other devices must be monitored to ensure that they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or a decrease in the device's electromagnetic immunity and lead to faulty operation.
- Failure to comply with the above can impair the performance characteristics of the device.

Initial use
The power supply unit is ready for use once it is plugged into a mains socket. The power supply unit must not be used outside of the specified input values as this may lead to the device being destroyed.

Technical data	
Type:	LXCPI2-006060BEH
Item number:	925 391
Input voltage:	100 – 240V AC, 50/60Hz
Max input current:	500mA
Output voltage:	6V DC, 600mA, only in combination with HARTMANN Veroval® blood pressure monitors
Max output current:	600mA
 Polarity of the DC voltage connection	
Output current protection:	Protection against short circuit
No-load losses:	< 0.075 watts
Approvals:	Medical standard IEC 60601-1 (in compliance with CISPR 11, IEC 61000-3-2, IEC 61000-3-3)
Conformity:	CE compliant as per European guidelines
Operating conditions:	Ambient temperature: 10 °C to 40 °C Relative humidity: 20 – 90 %, non-condensing
Storage / transport conditions:	Ambient temperature: -25 °C to +70 °C Relative humidity: 20 – 90 %, non-condensing
Dimensions L x B x H in mm:	62 x 39 x 26 (without mains plug)

Signs and symbols:	
	Caution
	Follow instructions for use
	Medical Device
	Do not use outdoors (indoor use only)
	Protective insulation/protection class 2
	Symbol for the marking of electrical and electronic equipment
	Cardboard recycling code

Subject to errors and changes

Veroval®

Keep dry

	Keep dry
	Batch code
	Catalogue number
	CE-marking according to European medical-device legislation
	Date of manufacture
	Manufacturer
	Authorised representative in the European Community (EC)
	Distributor
	Dispose of packaging in an environmentally responsible way
	Dispose of packaging in an environmentally responsible way
	Dispose of packaging in an environmentally responsible way
	Temperature limit
	Humidity limitation

Warranty conditions
■ We grant a 2-year warranty from the date of purchase for this high-quality power supply unit according to the conditions defined below.
■ Warranty claims must be asserted within the warranty period. The date of purchase must be proven by the purchase receipt.
■ Within the warranty period, HARTMANN will provide a replacement free of charge for any material and manufacturing defects on the device or will repair these defects. This does not result in an extension of the warranty period.
■ The device is only intended for the purpose described in these instructions for use.
■ Damage that has been caused by improper use or unauthorised operations are not covered by the warranty. Accessories that are subject to wear and tear are excluded from the warranty. Claims for damages are limited to the goods value; compensation for consequential loss is expressly excluded.
■ In the case of a warranty claim, please send the power supply unit with the purchase receipt to the responsible customer services department in your country either directly or via your dealer.

For a patient/user/third party in the European Union and in countries with identical regulatory regime: if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Packaging contents
1 Power supply unit
1 Operating instructions

Information on disposal of electronic devices
■ This product is subject to the European Directive 2012/19/EU on Used Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please obtain information about local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

Contact details for customer queries
ZA HARTMANN South Africa
Epsom Avenue, Northridge,
2194 Johannesburg
www.veroval.info
Tel. +27 860 4278 6266
ptzahelpdesk@hartmann.info

Made under the control of PAUL HARTMANN AG by
Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road
European Industrial Zone, Xiaolan Town
528415 Zhongshan City, Guangdong Province
P.R.C.

EC REP Donawa Lifescience Consulting Srl
Piazza Albania, 10
00153 Rome
Italy

Date of revision of the text: 2020-04-14 Version 3

	Caution
	Follow instructions for use
	Medical Device
	Do not use outdoors (indoor use only)
	Protective insulation/protection class 2
	Symbol for the marking of electrical and electronic equipment
	Cardboard recycling code

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Made under the control of PAUL HARTMANN AG by
Globacare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road
European Industrial Zone, Xiaolan Town
528415 Zhongshan City, Guangdong Province
P.R.C.

	Donawa Lifescience Consulting Srl <p>Piazza Albania, 10 00153 Rome Italy</p>
---	--

	Donawa Lifescience Consulting Srl <p>Piazza Albania, 10 00153 Rome Italy</p>
---	--

	Donawa Lifescience Consulting Srl <p>Piazza Albania, 10 00153 Rome Italy</p>
---	--

	Donawa Lifescience Consulting Srl <p>Piazza Albania, 10 00153 Rome Italy</p>
---	--

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