

SURFLO™-W


Winged I.V. catheter - FEP

SURFLO™ S-PLUS

Straight I.V. catheter - FEP

	<p>EN – Medical device IT- Dispositivo medico ES – Dispositivo médico FR – Dispositif Médical CS - Zdravotnický prostředek DA – Medicinsk udstyr DE – Medizinprodukt EL – Ιατροτεχνολογικό προϊόν ET – meditsiiniseade HU – Orvostechnikai eszköz LT - Medicinos prietaisai NO – Medisinsk utstyr NL – Medisch hulpmiddel</p>	<p>PL - Urządzenie medyczne PT – Dispositivo médico RO - DM TR – Tibbi cihaz SL - Medicinski pripomoček SV – Medicinteknisk produkt FI – Lääkinnällinen laite SK –Zdravotnícka pomôcka LV - Medicīniska ierīce BG - Медицински продукт RU – Медицинское устройство SR - Medicinski uređaj HR - Medicinski proizvod UK - Медичний пристрій</p>
	<p>EN – Distributor IT- Distributore ES – DISTRIBUIDOR FR – Distributeur CS - Distributor DA – Distributør DE – Händler EL – Διανομέας ET – edasimüüja HU – Forgalmazó LT - Platintojas NO – Forhandler NL – Distributeur</p>	<p>PL - Dystrybutor PT – Distribuidor RO - distribuitor TR – Distribütör SL - Distributer SV – Distributör FI – Jakelija SK –Distribútor LV - Izplatītājs BG - Дистрибутор RU – Дистрибьютор SR - Distributer HR - Distributer UK - Дистриб'ютор</p>
	<p>EN - Single sterile barrier system IT – Sistema barriera sterile singola ES – Sistema de barrera estéril individual FR – Système barrière stérile simple CS - Sterilní jednobariérový systém DA – Enkelt sterilt barrieresystem DE – Einfaches Sterilbarrieresystem EL – Σύστημα μονού φραγμού αποστείρωσης ET – Steriilne barjäärisüsteem HU – – Egyetlen steril zárórendszer LT - Vienkartinė sterili barjerinė sistema NO – Enkelt sterilt barrieresystem NL – Enkel steriel barrièresysteem PL - Pojedynczy sterylony system barierowy</p>	<p>PT – Sistema de barreira estéril única RO - Sistem de barieră sterilă unică TR – Tekli steril bariyer sistemi SL - Enojni sterilni pregradni sistem SV – Enkelt sterilt barriärsystem FI – yksi steriili estojärjestelmä SK – Systém jedinej sterilnej bariéry LV - Viena sterila barjeru sistēma BG - Единична стерилна система на бариера RU – Одноразовая стерильная барьерная система SR - Jednostruki sistem sterilne barijere HR - Jednostruki sustav sterilne barijere UK - Стерильна захисна система одноразового використання</p>

	<p>EN – Caution IT – Attenzione ES – Atención FR – Attention CS - Upozornění DA – Vigtigt DE – Achtung EL – Προσοχή ET – Tähelepanu HU – Figyelmeztetés LT - Dėmesio NO – OBS! NL – Opgelet</p>	<p>PL - Uwaga PT – Atenção RO - Atenție TR – Dikkat SL - Pozor SV – Obs! FI – Varoitus SK – Upozornenie LV - Uzmanību BG - Внимание RU – Внимание SR - Опрезност HR - Opresz UK - Увага</p>
	<p>EN – Do not reuse IT- Monouso ES – No reutilizer FR – Usage unique CS - K jednorázovému použití DA – Til engangsbrug DE – Einweg EL – Προϊόν μιας χρήσης ET – Ühekordne HU – Egyszer használatos LT - Vienkartinis NO – Til engangsbruk NL – Eenmalig gebruik PL - Do jednorazowego użytku</p>	<p>PT – Uso único RO - De unică folosință TR – Tek kullanımlık SL - za enkratno uporabo SV – Engångsbruk FI – Kertakäyttöinen SK – Na jednorazové použitie LV - Vienreizējai lietošanai BG - за еднократна употреба RU – Не использовать повторно SR - Za jednokratnu upotrebu HR - Nemojte ponovno upotrebljavati / samo za jednokratnu upotrebu UK- Одноразового використання</p>
	<p>EN – Do not re-sterilize ES – No re-esterilizar IT- Non risterilizzare FR – Ne pas re-stériliser CS – Neprovádějte opětovnou sterilizaci DA - Kan ikke gensteriliseres DE – Nicht erneut sterilisieren EL - Μην επαναποστειρώνετε ET – Mitte uuesti steriliseerida HU – Ne sterilizálja újra LT - Pakartotinai nesterilizuoti NO – Må ikke reesteriliseres NL - Niet hersteriliseren</p>	<p>PL - Nie sterylizować ponownie PT – Não reesterilizar RO - Nu reesterilizați TR – Yeniden sterilize etmeyin SL - Ne sterilizirajte ponovno SV - Återsterilisera inte FI – Ei saa steriloida uudelleen SK – Nevykonávajte opätovnú sterilizáciu LV - Nesterilizēt atkārtoti BG - Не стерилизирайте повторно RU – Не стерилизовать повторно SR – Nemojte sterilizovati HR - Nemojte ponovno sterilizirati UK - Не стерилізувати повторно</p>
	<p>EN – Consult instructions for use IT- Consultare le istruzioni per l'uso ES – Consulte instrucciones de uso FR – Suivez les instructions d'utilisation CS - Přečtěte si návod k použití DA - Se brugsanvisningen DE – Gebrauchsanweisung beachten EL - Διαβάστε τις οδηγίες χρήσης ET – Lugege kasutusjuhend läbi HU – Olvassa el a használati útmutatót LT - Žiūrėti naudojimo instrukcijas NO – Det henvises til bruksanvisningen NL - Raadpleeg de gebruiksaanwijzing</p>	<p>PL - Zapoznać się z instrukcją użytkowania PT – Consultar as instruções de uso RO - Citiți instrucțiunile de utilizare TR - Kullanma talimatlarına başvurun SL - Preberite navodila za uporabo SV - Läs bruksanvisningen FI – Katso käyttöohjetta SK – Prečítajte si návod na použitie LV - Saktīt lietošanas instrukcijas BG - Консултирайте инструкциите за употреба RU - См. инструкции по использованию SR - Pogledajte uputstva za upotrebu HR - Pročitajte upute za uporabu UK - Прочитайте інструкцію з використання</p>

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EN - INTENDED USE

Peripheral IV Catheter for short term peripheral venous access that, in combination with other medical device, allows the administration of fluids. Once placed in the vein, the IV Catheter can be connected to others medical devices to administer therapeutical solutions or drugs. IV Catheters can be connected with others medical devices provided with luer lock or luer fitting connections such as infusion sets, extension lines, stopcocks and syringes. In case of high pressure administration of fluids, like contrast media, the device is connected to extension lines connected to power injectors to administer the contrast media at high pressure.

The Catheter can be used on any patient population while taking into account the vascular anatomy of the patient and of the adequacy of the procedure.

The catheters are suitable for use with pressure injectors (max. 325psi)

DO NOT USE 26G IV CATHETERS FOR HIGH PRESSURE TREATMENTS.

CONTRAINDICATION

Do not use on patients with known hypersensitivity to the materials used.

MATERIALS USED

Stainless Steel, radiopaque Teflon (FEP), Polypropylene (PP), Polyoxymethylene (POM). Lubricants: polymethylsiloxane and aminofunctional polymethylsiloxane

Non pyrogenic

GENERAL WARNINGS

- Choose the correct IV Cannula size.
- Open the primary packaging following the indications (PEEL).
- Intended for single patient use only.
- Do not resterilize.
- DO NOT RE-USE. Re-use may lead to infection or other illness/injury and may compromise device functionality.
- Exposure to blood, either through percutaneous puncture with a contaminated needle or via mucous membranes or nonintact skin, may lead to serious illness such as hepatitis, HIV (AIDS), or other infectious diseases.
- Do not leave Luer slip connections unattended as they may become disconnected and result in blood exposure or blood loss.
- After use, discard the device according to your facility policy.
- Use the device to administer proper solutions and fluids in order to reduce as far as possible phlebitis caused by hypotonic or hypertonic solution.

GENERAL PRECAUTIONS

- For proper use, clinicians should be trained in vascular access and use of these catheters.
- Use these catheters in conjunction with your facility policy.
- Ensure aseptic technique, proper skin preparation, and continued protection of the site consistent with accepted standards of practice and your facility policy.
- Replace catheter according to your facility policy, relevant guidelines, or if the integrity of the device has been compromised.
- Follow Universal/Standard Precautions when inserting, maintaining, and removing these catheters to prevent exposure to bloodborne pathogens.
- Adhere to all the contraindications, warnings, precautions, and instruction for all infusate, as specified by their manufacturer
- Do not bend the needle before or during use.
- Do not use scissors or other sharp instruments at or near the insertion site.
- Use only with ISO 80369-7 compliant Luer connections. Non-ISO compliant Luer connections may cause leakage.
- Do not overtighten ISO 80369-7 Luer connections as damage may occur.

INSTRUCTIONS

1. **CAUTION:** Do not use if package is damaged or opened or if the expiration date has passed.
2. Prepare insertion site according to your facility policy.
3. Open package and remove catheter by grasping the grip.
4. Holding the wings (for winged IV catheter) or the catheter hub (for straight IV catheter) twist and remove needle cover in a straight, outward motion. Discard needle cover according to your facility policy.
5. Inspect device for damage before insertion. **CAUTION:** Do not use if device is damaged.
6. Access vessel using a low angle of insertion. **WARNING:** If the needle is partially or completely withdrawn from the catheter tubing during insertion, do not re-insert the needle into the catheter tubing as damage may occur.


7. Observe blood return in the flash chamber
8. Advance the catheter off the needle into the vessel. **CAUTION:** Do not pull back on the needle during catheter advancement.
9. Stabilize the catheter by gently pressing the wings to the skin (for winged IV catheter) or the catheter hub (for straight IV catheter) and apply venous compression. **WARNING:** Venous compression is required to prevent blood leakage from the catheter hub.
10. Pull the needle out of the catheter. **WARNING:** Immediately dispose of the needle in a puncture resistant, leak-proof sharps container keeping the needle point away from the body and finger at all times.
11. Continue to apply venous compression.
12. Remove Luer lock cap before discarding the needle into a puncture resistant, leak-proof sharps container.
13. Securely connect the Luer lock cap or other Luer device to the catheter hub.
14. Stabilize catheter and apply a sterile dressing according to your facility policy.
15. Flush or begin infusion.
16. Perform vascular access site care and dressing changes according to your facility policy.
17. Upon removal, examine catheter to ensure it is intact and discard according to your facility policy

Warnings (PRECAUTIONS.)


- Read the instructions before use.
- Use protective gloves.
- The product must be used immediately after the packaging has been opened.
- In case of incorrect transport and/or manipulation, the device or packaging could be subject to structural and/or functional damage.
- Do not use scissors at or near the insertion site.
- If the device is used at a high pressure or with injectors:
 - Connect directly to the pressure infusion system with the end luer-lock connector of the device.
 - Remove all the accessories connected to the device and replace them with a luer-lock cap where necessary.
 - Always check the patency of the device before use.
 - Never exceed the maximum pressure of 325 psi
- Flush the device immediately after the administration of medicines or biological fluids.
- Immediately remove any needle that has no coating, always keeping the tip away from your body and fingers.
- Do not expose to heat or direct sunlight.

Note: Any serious incident that has occurred in relation to the device should be reported to the Manufacturer and the Competent Authority of the Member State in which the user and /or the patient is established.

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